

MODUMATTRESS

**MODU**  
**FORM**®



*Shown with Moxie Platform Bed – MX9*

# MODUMATTRESS – VINYL TICKING

**MTF3680VMX-DD**

**36”d x 80”w x 6”h**

- Vinyl Ticking – non-fading and non-bleeding.
- Wipes clean with soap and water, or properly diluted disinfectant (final rinse with water to remove all residues).
- Cover treated with a fungistatic agent to protect from fungal growth and against degradation from micro-organisms and/or bacteria.
- Passes: AATCC Method 30-1988, AATCC Method 147-1998, AATCC Method 100-1999
- Skin irritation test Draize dermal result, non-allergenic.
- Meets Cal TB129.



- Bonded seams allows the mattress to bend without cracking.

- The dual-density foam core provides a greater range of comfort and support depending on the needs of the patients.

*Shown with Moxie Platform Bed – MX9*



# MODUMATTRESS – VINYL TICKING – INTEGRATED PILLOW



**MTF3680VMX-DD1P**  
1 Pillow



**MTF3680VMX-DD2P**  
2 Pillows



*Shown with Moxie Platform Bed – MX9*

# MODUMATTRESS – POLYURETHANE TICKING

**MTF3680MX-DD**  
**36”d x 80”w x 6”h**

- Thermal polyurethane coated nylon fabric - will not crack and is pinhole free.
- Wipes clean with soap and water, or properly diluted disinfectant (final rinse with water to remove all residues).
- Cover treated with a fungistatic agent to protect from fungal growth and against degradation from micro-organisms and/or bacteria.
- Passes: AATCC Method 30-1988, AATCC Method 147-1998, AATCC Method 100-1999
- Skin irritation test Draize dermal result, non-allergenic.
- Meets Cal TB129.



- Bonded seams allows the mattress to bend without cracking.

- The dual-density foam core provides a greater range of comfort and support depending on the needs of the patients.

*Shown with Moxie Platform Bed – MX9*



MODUMATTRESS – POLYURETHANE TICKING – INTEGRATED PILLOW



**MTF3680MX-DD1P**  
1 Pillow



**MTF3680MX-DD2P**  
2 Pillows



*Shown with Moxie Platform Bed – MX9*

# MODUMATTRESS – RIPSTOP TICKING

**MTF3680SAMX-DD**

**36”d x 80”w x 6”h**

- Ripstop fabric – Puncture resistant and rip-proof.
- Wipes clean with soap and water, or properly diluted disinfectant (final rinse with water to remove all residues).
- Cover treated with a fungistatic agent to protect from fungal growth and against degradation from micro-organisms and/or bacteria.
- Passes: AATCC Method 30-1988, AATCC Method 147-1998, AATCC Method 100-1999
- Skin irritation test Draize dermal result, non-allergenic.
- Meets Cal TB129.



- Bonded seams allows the mattress to bend without cracking.

- The dual-density foam core provides a greater range of comfort and support depending on the needs of the patients.

*Shown with Moxie Platform Bed – MX9*



# MODUMATTRESS – RIPSTOP TICKING – INTEGRATED PILLOW



**MTF3680SAMX-DD1P**  
1 Pillow



**MTF3680SAMX-DD2P**  
2 Pillows



*Shown with Moxie Platform Bed – MX9*

# MODU FORM<sup>®</sup>



## MODUMATTRESS

ModuForm, Inc.  
172 Industrial Road  
Fitchburg, MA 01420  
(978) 345-7942  
[moduform.com](http://moduform.com)