



FEATURES & BENEFITS

One-Piece Rotomolded Polyethylene Base & Arms

Seamless, monolithic construction eliminates joints, gaps, and weak points—preventing contraband concealment while resisting impact, fluids, and aggressive cleaning protocols.

Passive, Weight-Activated Tilt

Recline is engaged naturally through body weight and rearward pressure on the arms, eliminating levers, cords, or controls that could introduce ligature risk or require staff intervention.

Heavy-Duty Through-Bolted Recliner Mechanism

Steel recliner mechanism is mechanically fastened directly through the base, delivering exceptional in-use security, structural integrity, and long-term durability in high-use behavioral health environments.

Secure Yet Serviceable Mechanism Design

While rigidly fixed during use, the mechanism is intentionally engineered for controlled access and removal from the underside, allowing straightforward service or replacement without compromising safety or performance.

Integrated Clean-Out Zones

Open geometry around the base and arm structure allows full access for sanitation, supporting infection control protocols and minimizing fluid entrapment.

Upholstered Seat & Back for Supportive Comfort

Soft, supportive upholstery provides comfort appropriate for extended seated use while maintaining durability suited for continuously monitored care environments.

Stable, Low-Profile Base Geometry

Wide stance and grounded mass enhance stability during entry, exit, and recline, reducing tipping risk without relying on excessive weight or anchoring.

Designed for Continuously Monitored, Lower-Acuity Settings

Purpose-built for environments such as EMPATH units and crisis stabilization spaces where comfort, observation, and safety must coexist without over-engineering.

Observation Support Recliner

The Recliner Collection

SPECIFICATIONS

M100
43"d x 36"w x 46"h
19" Seat Height
28" Arm Height






Meets CAL TB 117-2013 flammability requirements
Meets MAS Green low-emissions criteria

**MODU
FORM**

172 Industrial Rd., Fitchburg, MA 01420
978.345.7942 | www.moduform.com

Observation Support Recliner

The Recliner Collection



**MODU
FORM**

172 Industrial Rd., Fitchburg, MA 01420
978.345.7942 | www.moduform.com